

Continuous Chemical Synthesis within the ResonantAcoustic® Mixing Platform

Joe Mayne, PhD
New Product and Process Development Manager
Resodyn Corporation



Technical
InterChange 2025

1

Objectives

- Understand the challenges in reactor design
- Describe behavior of liquid and gas mixing in a ResonantAcoustic® Continuous Flow Device
- Demonstrate benefits of using ResonantAcoustic® Continuous Flow Device for Chemical Reactions
 - Improved Mass Transfer (Decreased Mixing Time)
 - Increased Reaction Rate for Mixing Limited Reactions
 - Improved Heat Transfer
 - Able to process multiphase processes
 - Scalable Process



2

Technical
InterChange 2025

2

What Makes a Good Reactor Design?

- High Conversion, High Yield
- Fast, Complete Mixing of Reactants
- Temperature Control
- Consistent Results
- Flexible to range of operations
- All physical phenomena are well understood
- Physical phenomena are scalable

Example Rate Expression

$$r_P = k_1 C_A C_B$$

Optimized
Temperature

Optimized
Mixing




3

Technical
InterChange 2025

3


Legacy Batch vs. Advanced Flow Reactors

Chemglass 10L Reactor



Conventional Batch Reactors:
 + Inexpensive and Flexible
 + Multiphase Processes
 - Variability on product quality
 - Mixing and Heat-Exchange Limitations

Corning AdvancedFlow™ Reactor



Microfluidic Reactors:
 + Improved Mixing and Temp Control
 + Control and Tunability of Product
 + Decreased Waste
 - Complexity and Cost
 - Difficulty Processing Solids

4

Technical InterChange 2025

4

ResonantAcoustic® Mixing Phenomenon



11 inches

14 inches

Dow 200 @ 1,000 cP
90 g of acceleration


5

Technical InterChange 2025

5

Basic Principles of Gas-Liquid Mixing

- Faraday Instabilities Provide a Mechanism for Gas Bubble Encapsulation and Liquid Droplet Formation
- Once Droplets and/or Bubbles Are Created, Acoustic Waves Continue to Create Differential Motion Between the Two Phases
- Disparate Densities of Liquids and Gases Result in Excellent Mixing from RAM Motion




6

Technical InterChange 2025

6

Continuous Gas-Liquid Mixing Regimes

Compressive Regime vs. Splitting and Combining Regime



0 g Acceleration, 150 mL/min
4 mm Inner Diameter
Re = 780 → laminar flow

40 g Acceleration,
150 mL/min

100 g Acceleration,
150 mL/min

Gas Headspace is Required for RAM
Enhanced Mixing

Resodyn
ACOUSTIC MIXERS


7

Technical
InterChange 2025

7

Splitting and Combining Regime

Increases Surface Area and Reduces Boundary Layer Thickness




Resodyn
ACOUSTIC MIXERS

8

Technical
InterChange 2025

8

Demonstration of Minimal Back-Mixing



1% Xanthan Gum Solution (200 cP)
95g, 50 mL/min, 60 frames per second

Resodyn
ACOUSTIC MIXERS

9

Technical
InterChange 2025

9

Mass Transfer Coefficient Measurement

- Dissolved Oxygen Experiments:
 - Air-saturated water fed to 4mm reactor
 - Pure N₂ also fed at reactor inlet
 - Dissolved-oxygen meter measures transfer of O₂ into the gas phase

$$k_L a = \frac{1}{\tau} \ln \frac{c_{in} - c_{equ}}{c_{out} - c_{equ}}$$

10

10

Mass Transfer Coefficient Measurement

- Dissolved Oxygen Experiments:
 - Corning AFR is 7-fold improvement in Mass Transfer compared to legacy CSTR
 - Continuous ResonantAcoustic® Mixing is a further 35% improvement in Mass Transfer with an 88% lower residence time

Acceleration (g)	Configuration	k _L a (1/sec)
0	Corning AFR, 30% Gas	1.4
0	Corning AFR, 15% Gas	0.4
0	CSTR	0.1
60	RAM Flow-cell, 33% Gas	1.5
100	RAM Flow-cell, 10% Gas	0.9

11

11

Precipitation Reaction Demonstration

- Advanced microreactors easily clog in the presence of solids
- Continuous ResonantAcoustic® process can handle solids
- Example: Formation of Basic Copper Carbonate solids

$$2 Cu^{2+}(aq) + 2 CO_3^{2-}(aq) + H_2O \rightarrow Cu_2(OH)_2CO_3(s) + CO_2(g)$$

12

12

Precipitation Reaction Demonstration

Sample collected from microreactor 80-100 g acceleration

After Centrifuging

Solids Recovered from Centrifuge Vial

Microscope Images Reveal Basic Copper Carbonate Particles

13

Resodyn Acoustic Mixers InterChange 2025

13

Model Reaction Study: Mixing Time Measurement

- Parallel Competitive Reactions to quantify mixing time

$$H_2BO_3^- + H^+ \xrightarrow{\text{fast}} H_3BO_3$$

$$IO_3^- + 5I^- + 6H^+ \xrightarrow{\text{slower}} 3I_2 + 3H_2O$$

$$I_2 + I^- \leftrightarrow I_3^-$$

- Fast reaction consumes acid when perfectly mixed
- Any degree of poor mixing, side reaction will form I₂ formation
- UV/VIS measurement of [I₃⁻] in product, used to calculate mix quality and mixing time

14

Resodyn Acoustic Mixers InterChange 2025

14

Model Reaction Study: Mixing Time Measurement

Acceleration (g)	Corning AFR Mix Quality (%)	RAM Flow Cell Mix Quality (%)	Corning AFR Mixing Time (ms)	RAM Flow Cell Mixing Time (ms)
0	~65	~65	~65	~10
20	~75	~75	~45	~10
40	~85	~85	~25	~10
60	~90	~90	~15	~10
80	~92	~92	~10	~10
100	~95	~95	~5	~10

15

Resodyn Acoustic Mixers InterChange 2025

15

ResonantAcoustic® Continuous Microreactor

- RACMR Prototype:
 - Incorporation of up to three reactants
 - Immediate and complete mixing of reactants
 - RAM results in plug flow fluid profile
 - Jacketed for cooling or heating of reaction mixture ($4,400 \text{ W/m}^2\text{-K}$)

RACMR Design

16

Resodyn Acoustic Mixers

Technical InterChange 2025

16

Summary

- Rapid and complete incorporation of feeds (mixing time: 6 ms)
- Highly effective heat transfer ($4,400 \text{ W/m}^2\text{-K}$)
- Mass/Heat Transfer independent of flow velocity
- Ability to handle slurries, viscous flow, multi-phase processes, etc.
- Scalable process from bench-scale to industrial-scale reactors

17

Resodyn Acoustic Mixers

Technical InterChange 2025

17

Thank you for your time and attention.

18

Resodyn Acoustic Mixers

Technical InterChange 2025

18
